

ADVANCED CHEMICAL STORAGE & SAFETY CABINETS

Activated Charcoal Box

Filtration module adsorbing solvent vapour from cabinet extract air.

DESCRIPTION

A filter module that adsorbs solvent vapour from cabinet extract air, used where ducting to atmosphere is impractical. Filter life depends on what has been loaded onto it, so a change interval based on use rather than the calendar is what keeps it effective.

APPLICATIONS

Recirculating storage cabinets, ductless installations, retrofits in buildings that cannot be penetrated.

METHOD & COMPLIANCE REFERENCE

Filter saturation monitored and the element changed on a defined basis.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.