

PHARMACEUTICAL STANDARDS & SPECIALTY CHEMICALS

Analytical & Impurity Standards

Characterised impurity and degradation-product standards for method development.

DESCRIPTION

Characterised impurity, degradation product and related-substance standards, supplied with structural confirmation and assigned purity. Without them a related-substances method cannot be validated for specificity, and an unknown peak stays unknown.

APPLICATIONS

Related substances methods, forced degradation and stability-indicating method development, specificity validation, genotoxic impurity work.

METHOD & COMPLIANCE REFERENCE

Impurity control per ICH Q3A and Q3B; mutagenic impurities per ICH M7.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.