

PHARMACEUTICAL TESTING & ANALYTICAL SOLUTIONS

Bulk Density Tester

Measures the untapped bulk density of powders and granules.

DESCRIPTION

Determines the poured, untapped density of a powder or granule — the mass occupying a known volume when the material settles without compaction. It is the baseline against which tapped density is compared.

APPLICATIONS

Raw material acceptance, hopper and bin sizing, blend consistency checks.

METHOD & COMPLIANCE REFERENCE

Bulk density is described in USP <616> and Ph. Eur. 2.9.34.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.