

PHARMACEUTICAL STANDARDS & SPECIALTY CHEMICALS

Certified Reference Materials

Materials supplied with a certificate stating assigned values and uncertainty.

DESCRIPTION

Materials supplied with a certificate stating an assigned property value, its uncertainty and the traceability chain behind it, produced under an accredited reference material scheme. Uncertainty is the part most often skipped over — and the part an assessor asks about.

APPLICATIONS

Method validation, calibration, proficiency testing, measurement uncertainty budgets.

METHOD & COMPLIANCE REFERENCE

Produced under ISO 17034; certificates support the traceability requirements of ISO/IEC 17025.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.