

ADVANCED CHEMICAL STORAGE & SAFETY CABINETS

Combined Storage Cabinet

One cabinet with separated compartments for incompatible chemical classes.

DESCRIPTION

A single cabinet with separated compartments for incompatible classes — typically flammables on one side, acids and bases on the other — for laboratories where floor area does not allow separate units.

APPLICATIONS

Small laboratories, satellite storage points, teaching laboratories.

METHOD & COMPLIANCE REFERENCE

Flammable compartment per EN 14470-1; corrosive compartments segregated per safety data sheet.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.