

LABORATORY INFRASTRUCTURE & TECHNICAL FURNITURE

Emergency Shower / Eye Wash

Immediate decontamination equipment sited within reach of the hazard.

DESCRIPTION

Drench shower and eyewash sited within ten seconds of unobstructed travel from the hazard, delivering tepid water at the required flow for a fifteen-minute flush. Weekly activation is what keeps the line clear and the record complete.

APPLICATIONS

Wherever corrosives or solvents are handled — wet chemistry benches, dispensing, sampling booths, chemical stores.

METHOD & COMPLIANCE REFERENCE

Emergency eyewash and shower equipment per ANSI/ISEA Z358.1.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.