

WASTE & SAFETY HANDLING SYSTEMS

EZwaste® Cap Assembly

Sealed cap assembly converting a standard waste container into a closed system.

DESCRIPTION

A sealed cap that turns a standard waste container into a closed system, with ports for instrument waste lines and a carbon filter on the vent. Solvent vapour stays in the bottle instead of escaping into the laboratory through an open funnel.

APPLICATIONS

HPLC and LC-MS waste collection, sample preparation waste, general solvent collection at the bench.

METHOD & COMPLIANCE REFERENCE

Supports closed-system solvent waste handling and laboratory vapour control practice.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.