

WASTE & SAFETY HANDLING SYSTEMS

EZwaste® Heavy Duty Solvent Waste System

Higher-capacity closed solvent waste handling for demanding installations.

DESCRIPTION

Higher-capacity closed collection for laboratories generating solvent waste faster than a standard bottle can absorb, with the same sealed cap, filtered vent and level indication at larger volume.

APPLICATIONS

High-throughput chromatography suites, preparative HPLC, contract laboratories.

METHOD & COMPLIANCE REFERENCE

Closed-system handling at larger volume; secondary containment sized to the container.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.