

WASTE & SAFETY HANDLING SYSTEMS

EZwaste® Manifold

Multi-port manifold collecting several instrument waste lines into one container.

DESCRIPTION

Collects several instrument waste lines into a single container through one sealed manifold, so a bank of instruments shares one bottle rather than each carrying its own.

APPLICATIONS

Multi-instrument HPLC benches, shared preparation areas, automated systems with several waste streams.

METHOD & COMPLIANCE REFERENCE

Closed-system collection; containers filled only to the labelled working level.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.