

WASTE & SAFETY HANDLING SYSTEMS

EZwaste® UN/DOT Solvent Waste System

Solvent waste collection in UN/DOT-rated containers for transport compliance.

DESCRIPTION

Collection built on UN/DOT-rated containers, so the container the waste is collected in is also the container it can legally be transported in. That removes the decanting step, which is where most spills and exposures actually happen.

APPLICATIONS

Sites where solvent waste is collected by a licensed contractor, laboratories under transport-regulated disposal.

METHOD & COMPLIANCE REFERENCE

Containers UN/DOT rated for transport of the waste stream collected.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.