

**ADVANCED CHEMICAL STORAGE & SAFETY CABINETS**

# Flammable Safety Cabinet

Fire-rated cabinet for storing flammable liquids within the laboratory.

**DESCRIPTION**

Fire-rated storage that keeps flammable liquids usable inside the laboratory while giving the fire service a defined period before the contents become involved. Ratings are stated in minutes — 15, 30, 60 or 90 — and the rating you need follows from the building fire strategy, not from the volume stored.

**APPLICATIONS**

Solvent storage at the point of use, under-bench and freestanding installations, flammable stores.

**METHOD & COMPLIANCE REFERENCE**

Safety storage cabinets for flammable liquids per EN 14470-1, classified Type 30, 60 or 90.

**SUPPLIED AND SUPPORTED IN BANGLADESH**

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.