

PHARMACEUTICAL TESTING & ANALYTICAL SOLUTIONS

Laboratory Glasswashers & Dryers

Principal: Miele

Validated, reproducible cleaning and drying of laboratory glassware.

DESCRIPTION

Washes, rinses and dries laboratory glassware to a repeatable cycle, with injector racks that reach inside narrow-necked vessels. A conductivity-monitored final rinse and a captured cycle record are what make the cleaning defensible rather than merely done.

APPLICATIONS

QC and microbiology glassware reprocessing, volumetric glassware, cleaning validation support.

METHOD & COMPLIANCE REFERENCE

Cycles validated against a residue acceptance criterion; final rinse water quality per the applicable water monograph.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.