

PHARMACEUTICAL TESTING & ANALYTICAL SOLUTIONS

Leak Test Apparatus

Checks seal integrity of blister packs, strips and sealed containers.

DESCRIPTION

Places sealed packs under vacuum in a dye bath and looks for ingress, revealing openings the eye cannot find. Failures trace back to sealing temperature, pressure or tooling, which is what makes it useful as an in-process control rather than only a release test.

APPLICATIONS

Blister and strip seal integrity, in-process packaging control, ampoule and vial checks.

METHOD & COMPLIANCE REFERENCE

Container closure integrity approaches are discussed in USP <1207>; dye ingress is a probabilistic method within it.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.