

PHARMACEUTICAL STANDARDS & SPECIALTY CHEMICALS

Secondary & Working Standards

Standards qualified against a primary reference for routine analytical use.

DESCRIPTION

Standards qualified against a pharmacopoeial primary reference and supplied with the certificate showing that comparison. They carry the routine analytical load at a fraction of the cost and lead time of the primary standard, which is then reserved for qualifying the next working batch.

APPLICATIONS

Routine assay and content uniformity, system suitability, method transfer, day-to-day QC.

METHOD & COMPLIANCE REFERENCE

Traceable to USP, EP or BP primary reference standards; the qualification protocol accompanies the certificate.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.