

LABORATORY INSTRUMENTATION & PROCESS SYSTEMS

Thermostatic Water Bath

Holds samples at a controlled temperature for incubation and reaction work.

DESCRIPTION

Holds a bath of water at a set temperature with circulation and PID control, so samples sit in a stable, uniform environment for as long as the method requires. Uniformity across the working volume, not the setpoint on the display, is what determines whether two positions give the same result.

APPLICATIONS

Sample incubation, dissolution medium pre-heating, reagent warming, viscosity and solubility work.

METHOD & COMPLIANCE REFERENCE

Temperature uniformity and stability are mapped and qualified as part of IQ/OQ.

SUPPLIED AND SUPPORTED IN BANGLADESH

- Import, customs, delivery, site preparation and installation.
- IQ/OQ qualification, scheduled calibration and validation support.
- Genuine spare parts, consumables and operator training in Dhaka.

Specifications, capacities and model options vary by configuration and are confirmed with each quotation. This sheet is generated from our product catalogue; the full manufacturer datasheet is available on request.